 **CERVICAL-Stim[®]**
Osteogenesis Stimulator

U.S.A. EDITION


Table of Contents

Cervical-Stim® Instruction Manual

Glossary	1
Prescription Information	3
Adverse Events	3
Device Description	5
Device Life	6
Treatment Time	6
Charging/Recharging the Battery	7
Device Operation	8
Wearing the Device	9
Care and Cleaning	10
Travel	10
Storage	10
Disposal	10
Service	10
Visual and Audio Indicators	11
Equipment Classification and Device Symbol Descriptions	12
Clinical Data Summary	13
Warranty Information	15

To learn more about Orthofix, please visit our website at www.orthofix.com.

Package Contents:

- I- Cervical-Stim Osteogenesis Stimulator
- I- Literature Pack
- I- Power Supply
- I- Line Cord
- I- Comfort Collar

**THIS DEVICE IS NONSTERILE.
IT DOES NOT REQUIRE STERILIZATION.**

Glossary

adjunct/adjunctive: an additional therapy soon after the surgery

allograft: bone graft taken from another human

CT scan: Computerized Tomography Scan; a picture of structures within the body created by a computer that takes the data from multiple x-ray images and turns them into pictures on a screen

defibrillator: an electrical device used to restore normal heartbeat by applying a brief electric shock

diabetes mellitus: a chronic disease in which the body has a lack of insulin and cannot make use of carbohydrates

electrocardiogram: a recording of the electrical activity of the heart

epilepsy: also known as “seizure disorder”; a pattern of seizures caused when the nerve cells in the brain fire electrical impulses at a rate of up to four times higher than normal

fusion: the joining of two bones into a single unit and thus eliminating motion between the two bones

high-risk patients: patients who have conditions or health concerns that may hinder normal bone healing/fusion; these can include smoking, bone depleting medications, osteoporosis, multi-level fusions, revision surgery, etc.

ligamentous spinal trauma: typically occurring in the cervical spine; generally brought about by extreme forces causing sudden forward or backward bending of the spine (e.g., whiplash)

metastatic cancer: type of cancer that spreads (i.e., metastasizes) from its original site to another area of the body

MRI: Magnetic Resonance Imaging; procedure that uses a magnet connected to a computer to create images of internal structures of the body

noninvasive: does not require placing an instrument or device through the skin or body orifice to get treatment

nonunion fracture: failure of the pieces of a broken bone to bond together

osseous: having to do with bone, consisting of bone, or resembling bone

osteoporosis: progressive disease that causes thinning of the bones with a decrease in bone mass

pacemaker: surgically implanted electronic device used to stimulate or regulate contractions of the heart muscle

Paget's disease: a chronic bone disorder that typically results in enlarged, deformed bones

PEMF: Pulsed ElectroMagnetic Field; a time-varying (pulsing) magnetic field; an external coil that produces a magnetic field through which an electrical charge is given

postoperative: following an operation

prospective randomized controlled clinical study: a forward-looking study involving human beings in which half of the patients are, at random, given a treatment and half are not

renal disease: permanent loss of most of the kidneys' ability to remove waste and maintain fluid and chemical balance in the body

rheumatoid arthritis: a chronic disease resulting in stiffness and inflammation of the joints, loss of mobility, weakness and deformity

spondylitis: inflammatory disease of the spine; inflammation of one or more vertebrae of the spine

Prescription Information

Indication

The Cervical-Stim® is a noninvasive, pulsed electromagnetic bone growth stimulator indicated as an adjunct to cervical fusion surgery in patients at high-risk for non-fusion.

Contraindications

There are no known contraindications for the Cervical-Stim as an adjunct to cervical spine fusion surgery.

Warnings

- Do not use Cervical-Stim if you have a cardiac pacemaker or defibrillator because it may interfere with the operation of your pacemaker or defibrillator. If you use the Cervical-Stim and it affects your pacemaker or defibrillator, it may injure your heart. Consult your cardiologist.
- Remove the Cervical-Stim prior to any imaging procedures (e.g., CT scan, MRI, etc.). If you wear the Cervical-Stim during these procedures, you could be injured, the imaging being produced may be ruined, and/or the Cervical-Stim could be damaged.

Precautions

- Avoid using the Cervical-Stim if you do not understand the instructions your doctor has given you. If you use the Cervical-Stim incorrectly, it may harm you or may not help your healing process.
- The Cervical-Stim has not been evaluated in treating patients with the following conditions: osseous or ligamentous spinal trauma, spondylitis, Paget's disease, moderate to severe osteoporosis, metastatic cancer, renal disease, rheumatoid arthritis, uncontrolled diabetes mellitus, patients prone to vascular migraine headache, seizure, epilepsy, thyroid conditions or neurological diseases.
- Animal reproductive studies performed with this device did not show any harmful effects in animals. However, the safety of this device for use on patients who are pregnant or nursing has not been established.

Adverse Effects

Adverse effects may be experienced when using the Cervical-Stim. These adverse effects may include: increased pain, numbness and tingling, headache, migraines and nausea. These effects may or may not be directly related to use of the Cervical-Stim. Any adverse effects that are related to the Cervical-Stim should stop when you discontinue use. The following table shows all adverse events that were reported during the clinical study.

Adverse Events

Adverse Events Reported at 6 Months by Treatment Group

Adverse Events	Control Group (n=160)		Cervical-Stim Group (n=163)	
	# (%) of Events	# (%) ¹ of Patients Experiencing the Event	#* (%) of Events	# (%) ¹ of Patients Experiencing the Event
Increased Neck Pain	10 (14.9)	9(5.6)	16(17.8)	15(9.2)
Shoulder/Arm Pain	10(14.9)	9(5.6)	16(17.8)	16(9.8)
Re-Injury to Cervical Spine	10(14.9)	8(5.0)	9(10.0)	9(5.5)
Adjacent level pathology	3(4.5)	3(1.9)	8(8.8)	8(4.9)
Surgical Complications	2(3.0)	2(1.3)	7(7.7)	5(3.1)
LBP/Lumbar pathology	8(11.9)	8(5.0)	5(5.5)	5(3.1)
Trauma/Injury(not cervical)	2(3.0)	2(1.3)	5(5.5)	4(2.5)
Numbness/Tingling	6(8.9)	6(3.8)	4(4.4)	4(2.5)
Headache/Migraine	2(3.0)	2(1.3)	4(4.4)	4(2.5)
Nonspecific/Unrelated Pain	2(3.0)	2(1.3)	3(3.3)	3(1.8)
Nausea	0	0	2(2.2)	2(1.2)
Dizziness/Vertigo	2(3.0)	2(1.3)	1(1.1)	1(0.6)
Rash/Discoloration	0	0	1(1.1)	1(0.6)
Rapid/Irregular Heartbeat	0	0	1(1.1)	1(0.6)
Shortness of Breath	0	0	1(1.1)	1(0.6)
Ringing in Ears	0	0	1(1.1)	1(0.6)
Neurologic Symptom/Stroke	1(1.5)	1(0.6)	1(1.1)	1(0.6)
Lump in Throat	0	0	1(1.1)	1(0.6)
Diagnosis of Diabetes	0	0	1(1.1)	1(0.6)
Diagnosis of Breast Cancer	0	0	1(1.1)	1(0.6)
Seizure	0	0	1(1.1)	1(0.6)
Death, Unrelated	0	0	1(1.1)	1(0.6)
Tenderness	1(1.5)	1(0.6)	0	0
Screw Broken	1(1.5)	1(0.6)	0	0
Graft Collapse	1(1.5)	1(0.6)	0	0
Carpal Tunnel Syndrome	2(3.0)	2(1.3)	0	0
Choking Sensation	1(1.5)	1(0.6)	0	0
Cardiac Symptoms	1(1.5)	1(0.6)	0	0
Nephrotic Syndrome	1(1.5)	1(0.6)	0	0
Suicide Attempt	1(1.5)	1(0.6)	0	0
TOTAL	67	47 ²	90	58 ²

¹ % expressed as number of patients experiencing the event / total number of patients in the group

² Some patients experienced multiple adverse events

* There were several adverse events that were more frequently observed in the Cervical-Stim group than in the control group. Given the types of events, it is unlikely that these adverse events are related to the treatment.

For more information regarding safety, refer to the section entitled “Clinical Data Summary.”

Device Description


The Cervical-Stim is an external, low-level, pulsed electromagnetic field device intended to be worn after your cervical fusion surgery to increase your chances of a successful fusion. The Cervical-Stim has been designed with your comfort and convenience in mind. It is a single-piece device that is lightweight, flexible and portable. It allows you to move freely during your treatment. A Liquid Crystal Display (LCD) and an audible alarm give you information during your treatment such as the operational status, treatment time remaining, battery capacity, etc. See “Visual and Audio Indicators” for more information.

The Cervical-Stim is comprised of a control unit and a treatment transducer. The control unit contains a computer chip that generates the electrical signal. That signal is converted to a uniform, low-energy magnetic field by the treatment transducer. By placing the device directly over the treatment area, the PEMF signal is delivered directly to your fusion site.

The Cervical-Stim is powered by a rechargeable lithium-ion battery pack. The LCD and audible alarm will alert the patient when the battery is low and needs to be recharged. See “Charging/Recharging the Battery” for more information.

To ensure that the device is functioning properly, the Cervical-Stim constantly monitors battery voltage and the electrical signal. If at any time during treatment the device stops functioning properly, the LCD will display an appropriate symbol or error code. See “Visual and Audio Indicators” for more information.

Cervical-Stim Osteogenesis Stimulator - Model 2505

Device Life

The Cervical-Stim can provide you with up to 270 daily treatments of four hours each. Your doctor will decide how long (weeks/months) your overall length of treatment will be.

Treatment Time

The Cervical-Stim should be worn for four hours each day. The device will automatically turn off when four hours of treatment is reached in a 24-hour day. The Cervical-Stim may be used at any time of day that is convenient. You may turn the device off at any time by simply pressing the On/Off button on the control panel. It is lightweight and adjustable. And because the Cervical-Stim is portable, you can get treatment while you are sitting, walking, reclining, sleeping, etc. However, since each patient is unique, your activity level should be based on your doctor's instructions. Not following your doctor's instructions can cause injury or diminish the effect of wearing the Cervical-Stim.

Charging/Recharging the Battery


The Cervical-Stim is powered by a rechargeable lithium-ion battery pack. A power supply is provided with the device. Use only the Orthofix charging system to charge the battery.


Note: The Cervical-Stim battery will require charging prior to the first use.

To charge/recharge the battery, simply plug the barrel connector end of the power supply into the charger port located on the control unit. Plug the line cord securely into the power supply. Plug the line cord into any standard AC wall outlet. A fully discharged battery may require up to 4 hours to charge completely.

The Cervical-Stim battery can be recharged at any time the device is not in use. It is strongly recommended that the device be recharged after completing daily treatment.

Note: The Cervical-Stim will not deliver treatment while charging.

When the device is on, the Cervical-Stim LCD will show a battery capacity symbol. A flashing battery outline, the symbol  and an audible beep indicate a battery low condition and that the battery requires recharging. See “Visual and Audio Indicators” for more information.

Device Operation

Turning the Device On and Off

The Cervical-Stim is turned on and off by pressing the On/Off button on the control unit of the device. When the device is on, a sequence of status messages will display momentarily. The LCD should then show the treatment time remaining and a flashing Orthofix logo. The flashing logo indicates that the device is on and functioning normally. (If you do not see this on the display, contact your provider.) A backlight button is on the control unit. In low light, press the backlight button for illumination of the LCD.

Timing of Treatment Sessions

To stop treatment prior to the end of the daily treatment session, simply press the On/Off button. To resume treatment, press the On/Off button again. The LCD will display the remaining treatment time.

Note: For the countdown to function correctly, treatment sessions should be greater than 60 minutes duration.

Wearing the Device

The Cervical-Stim is intended for the cervical spine and may be worn with or without a brace. For better comfort, clothing should be worn between the skin and the Cervical-Stim.

To use the Cervical-Stim, simply slip the device over the head so that it rests comfortably against the neck and shoulders (see figure below). Or, the device may be opened at the buckle and placed against the neck and shoulders. Fasten the buckle like a seatbelt, to secure the unit.

A Comfort Collar is included with Cervical-Stim for improved fit and support. The Comfort Collar can be inserted or removed from the inside collar of the device using Velcro Tabs, which attach directly to the material.

Care and Cleaning

The Cervical-Stim is a sophisticated electronic device and should be handled with care. Dropping or other mistreatment of the Cervical-Stim may cause damage to the device.

DO NOT expose the Cervical-Stim to direct sunlight for long periods of time. The device control unit may be damaged.

DO NOT expose the Cervical-Stim to excessive heat. In warm climates, the temperature in a car or trunk can exceed 160°F / 71°C. Excessive heat can damage the control unit of the device.

DO NOT expose the Cervical-Stim to excessive moisture. Moisture can damage the electronic components of the device and the device may stop working.

DO NOT dispose of the Cervical-Stim in an incinerator.

DO NOT use solvents to clean the Cervical-Stim. Clean the device by wiping with a soft, damp cloth.

Travel

When traveling by air, it is best to check the Cervical-Stim with the luggage. If the device is taken on board the airplane, it should not be worn when passing through passenger screening devices. The Cervical-Stim could be damaged. The Cervical-Stim user manual should be taken with you to quickly and easily identify the device for any security personnel.

Storage

The Cervical-Stim should be stored within 14°F to 113°F (-10° C to 45° C) The Cervical-Stim operating temperature range should be within 41°F to 104°F (+5° C to 40° C). Relative Humidity: Up to 95%, non-condensing

Disposal

The Cervical-Stim is for single patient use. Product contains lithium batteries; do not incinerate. Dispose of device properly to prevent injury. Please dispose of this product at collection facilities for waste electrical equipment used in household.

Service

If you have questions concerning the device or require any assistance, please call 800-535-4492 or 214-937-2000. There are no user serviceable parts. Notify the manufacturer for any servicing needs.



Visual and Audio Indicators

The LCD and audible alarms are designed to provide helpful information to the user. The chart below shows the various displays and alarms and their meaning.

Cervical-Stim LCD Visual and Audio Indicators

Symbol /Alarm	Description	Meaning
	All LCD symbols visible and continuous audible alarm for approximately 5 seconds	Power-on self test
	Countdown timer displays remaining treatment time (hours & minutes) Orthofix logo flashes	Normal treatment in progress
	Countdown timer displays three dashes audible alarm (5 beeps)	No treatment time remaining
	Steady symbol for approximately 5 seconds	Treatment complete — power off
	Symbols flash / audible alarm (approximately 1 beep per second)	Battery low — recharge required
	Steady symbol indicates approximate % of charge Symbol filling repeatedly indicates charge mode	Battery status — remaining charge or charging mode
	Continuous audible alarm	Device locked — call for service
	Display of any E code (e.g., E01, E02 . . .)	Error message — call for service

Equipment Classification and Device Symbol Descriptions

Symbol	Meaning
	Attention - Refer to Instructions for Use
	Type BF Applied Part
	On/Off
	Backlight Button
	Storage Temperature Range
	Year of Manufacture for Active Device
	Charger Port

Equipment Classifications

- Internally powered equipment
- Type BF applied part
- IEC 529 enclosure rating: IPX0
- Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or nitrous oxide
- Mode of operation: intermittent operation

The use of accessories other than those specified may result in increased emissions or decreased immunity of the device.

The battery charger is provided with a 3-wire appliance inlet but considered double insulated with Class II construction throughout.

For safe usage, follow manufacturer instructions when using the product.

Use of the product in any other manner could have harmful effects and/or void the warranty.

Note: Inspect the device prior to each use for wear or deterioration.
Do not use if the device does not appear to be in suitable condition.

Clinical Data Summary

The Cervical-Stim was studied in humans to evaluate its safety and effectiveness as a therapy added to routine care (adjunct therapy) for high-risk patients having a cervical fusion surgery for degenerative conditions. Patients were high-risk if they were a smoker (one pack per day or more) and/or had a multi-level fusion surgery (more than one level).

The 323 patients were randomly assigned, similar to a coin toss, to one of two groups: either the control group (routine care only) or the treatment group (Cervical-Stim + routine care). One hundred and sixty (160) patients were assigned to the control group and 163 patients were assigned to the Cervical-Stim group. Patients wore the Cervical-Stim unit for 4 hours each day either for 4 continuous hours or in one hour sessions.

Safety and effectiveness was evaluated by measuring the following:

- rate and severity of adverse events
- rate of cervical fusion by six months after surgery as determined by x-ray

Eighty-four percent (84%) of the Cervical-Stim group were fused by six months (102/122 patients) versus only 69% of the control group (81/118 patients). This is a 15% difference between these two groups and is statistically significant (meaningful); $p=0.0065$. That is, more patients fused in the Cervical-Stim group than in the control group.

The rate of patients who came back for their six month examinations and x-rays was 74% for the Cervical-Stim group and 73% for the control group. Patients who did not come back for scheduled examinations could not be evaluated; thus their success or failure is not known. These unavailable data could have a positive or negative effect on the overall success of this study.

One hundred and twelve (112) patients reported a total of 157 adverse (negative) effects for both groups combined at six months after surgery. There was no significant (meaningful) difference in the total number of adverse events or the number of patients reporting effects in the control group and the Cervical-Stim group nor in the numbers of patients in each group who experienced an adverse event. The adverse effects that may be experienced include: increased pain, numbness and tingling, headache, migraines and nausea. These effects may or may not be directly related to the use of the Cervical-Stim.

Clinical success with regard to symptoms was evaluated by the following:

- no worsening in neurological function
- an improvement in pain
- no worsening in Neck Disability Index

Based on the criteria above, there was no major difference between the control group and the Cervical-Stim group in clinical success. An equal number of patients in both groups showed an improvement in their clinical condition after surgery, regardless of treatment.

Long-term x-ray information collected at 11 months after surgery or later showed no meaningful difference in fusion rate between the Cervical-Stim treatment group and the control group who received routine care alone.

The results of this study show that the use of the Cervical-Stim is both safe and effective in increasing the frequency of fusion by six months after surgery in high-risk subjects having cervical fusion.

Warranty Information

Orthofix Inc. warrants the Cervical-Stim to be free from defects in materials and workmanship for one year from the date of first use. Provided that all terms and conditions of this Limited Warranty are complied with, Orthofix Inc. will replace defective components.

This Limited Warranty applies to the product only under normal use and does not cover any damage or defect caused by accident, misuse, abuse, fire, flood, and acts of God or by any alteration, tampering, repair or attempted repair by anyone other than Orthofix Inc. This warranty only applies to the patient for whom the product is prescribed and is not assignable or transferable. Defective products covered by this Limited Warranty must be returned to Orthofix Inc. Attention: Orthofix Returns. You must call a Patient Services Representative at 1-800-535-4492 or your local distributor to obtain the Return Authorization (RA) number and address prior to returning the product.

Except as specifically required by applicable law, the foregoing warranty is in lieu of all other warranties, expressed or implied and Orthofix Inc. specifically disclaims any and all warranties of merchantability or fitness for a particular purpose. Under no circumstances shall Orthofix Inc., its authorized representative, affiliated or subsidiary companies be liable for special, consequential or incidental damages. The sole remedy with respect to any defective product shall be limited to replacement.

This Limited Warranty may not be extended or modified except in writing by Orthofix Inc. No sales person, representative, distributor or doctor is authorized to make or consent to any extension or modification of the terms of this Limited Warranty.

For additional information and/or device assistance, contact Orthofix Patient Services at 800-535-4492.

Manufactured by:
Orthofix Inc.
3451 Plano Parkway
Lewisville, TX 75056-9453
214-937-2000

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CS-1404 PL-US © Orthofix Holdings, Inc.