PHYSIO-Stim[®]
Osteogenesis Stimulator

U.S.A. EDITION


Table of Contents

Physio-Stim® Instruction Manual

Prescription Information	1
Orthofix PEMF Stimulation	2
Clinical Success of the Physio-Stim	2
Device Description	3
Treatment Instructions	4
Device Operation	4
Turning the Device On and Off	4
Timing of Treatment Sessions	5
Charging/Recharging the Battery	5
Visual and Audio Indicators	6
Device Application	7
Models 3202 and 3303	7
Model 3313.....	8
Models 3314L and 3314R	9
Model 3315	10
Equipment Classification and Device Symbol Descriptions	11
Care and Cleaning	12
Travel	12
Storage	12
Disposal	12
Service	12
Warranty Policy.....	13

Package Contents:

- I- Physio-Stim Osteogenesis Stimulator
- I- Literature Pack
- I- Power Supply
- I- Line Cord

THIS DEVICE IS NONSTERILE.
IT DOES NOT REQUIRE STERILIZATION.

U.S. Patent No. 5,743,844

Prescription Information

Indication

The Physio-Stim[®] is indicated for the treatment of an established nonunion acquired secondary to trauma, excluding vertebrae and all flat bones, where the width of the nonunion defect is less than one-half the width of the bone to be treated. A nonunion is considered to be established when the fracture site shows no visibly progressive signs of healing.

Contraindication

Use of this device is contraindicated where the individual has synovial pseudarthrosis.

Warnings

- The safety and effectiveness of the use of this device on individuals lacking skeletal maturity has not been established.
- In the presence of a malaligned nonunion, careful consideration of the use of this device must be undertaken on an individual basis, as treatment with this device is not intended to alter or affect the degree of malalignment.
- Demand type pacemaker operation may be adversely affected by exposure to pulsed electromagnetic fields. Physicians should not prescribe a Physio-Stim for application which may place the treatment transducer in close proximity to the pacemaker. Further screening by the attending cardiologist is recommended (such as with an electrocardiogram).
- Animal studies conducted to date do not suggest any long-term adverse effects from the use of this device. However, long-term effects in humans are unknown.
- The safety and effectiveness of this device on individuals with a nonunion secondary to, or in connection with, a pathological condition has not been established.

Precautions

- Nonunion fractures with gaps in excess of 1 centimeter (cm) have not been evaluated.
- Although animal reproductive studies performed with this device demonstrated no adverse findings, the safety of use of this device during pregnancy and nursing in humans has not been established.
- This device should not be used if there are mental or physical conditions which preclude patient compliance with physician and device instructions.

Adverse Events

Rare instances of reversible minor discomfort have been reported. They were: cumbersome or uncomfortable, tingling or pain and minor skin rash.

Orthofix PEMF Stimulation

Pulsed electromagnetic field (PEMF) osteogenesis stimulation is a safe, nonsurgical, prescription treatment to heal nonunion fractures and promote spinal fusion. Electrical currents have been used to heal bones since the mid-1800s. However, it wasn't until the 1950s that scientists made an important discovery. When human bone is bent or broken, it generates an electrical field. This low-level electrical field activates the body's own repair mechanism which, in turn, stimulates bone healing.

Orthofix PEMF osteogenesis stimulators generate a uniform, low-level, pulsed electromagnetic field similar to the electrical field generated by the body. The application of PEMF directly to the fracture site helps activate and augment the body's natural healing process to enhance bone fusion. Thank you for including Orthofix in your healing process. To learn more about osteogenesis stimulation, please visit our website at www.bonestimulation.com.

Clinical Success of the Physio-Stim


The Physio-Stim has been proven safe and effective in clinical studies. In a prospective, multicenter, controlled clinical study for nonunion fractures, an overall success rate of 80% was demonstrated among 126 patients (135 nonunion fractures) who averaged 3 or greater hours of daily treatment. The average duration of nonunion in these difficult fractures was 2.6 years, with an average of two prior surgical procedures per fracture. The success rate of Physio-Stim treatment for nonunion repair demonstrated no statistically significant change over long-term (four year) follow-up.

Device Description

The Physio-Stim is an external, low-level, PEMF device and has been designed with patient comfort and convenience in mind. It is a single-piece device that is lightweight, flexible and portable, allowing freedom of movement during treatment. A Liquid Crystal Display (LCD) and audible alarm provide important feedback during treatment such as the operational status, treatment time remaining, battery capacity, etc. See “Visual and Audio Indicators” for more information.

The Physio-Stim is comprised of a control unit and treatment transducer. The control unit contains a micro-processor that generates the Physio-Stim electrical signal. That signal is converted to a highly uniform, low-energy magnetic field by the treatment transducer. When the device is centered over the treatment area, the therapeutic Physio-Stim PEMF signal is delivered directly to the treatment site.

Model: 3313

The Physio-Stim is powered by a rechargeable lithium-ion battery pack. The LCD and audible alarm will alert the patient when the battery is low and needs to be recharged. See “Charging/Recharging the Battery” for more information. To ensure that the device is functioning properly, the Physio-Stim constantly monitors battery voltage and the electrical signal. If at any time during treatment the device stops functioning properly, the LCD will display an appropriate symbol or error code. See “Visual and Audio Indicators” for more information.


Treatment Instructions

The Physio-Stim should be worn for three hours per day. Based upon clinical data for nonunion, the overall treatment duration ranges between 90 and 180 days based upon specific patient conditions. At the end of the daily treatment the device will turn itself off. The Physio-Stim may be used at any time of day that is most convenient and comfortable for the patient. It is lightweight and adjustable. And because the Physio-Stim is portable, treatment can be received while the patient is sitting, walking, reclining, sleeping, etc. However, since each patient is unique, the overall activity level should be based on physician instructions.

Device Operation

Turning the Device On and Off

The Physio-Stim is turned on and off by pressing the On/Off button on the control unit of the device. When the device is on, a sequence of status messages will display momentarily. The LCD should then show the treatment time remaining and a flashing Orthofix logo. The flashing logo indicates that the device is on and functioning normally. (If you do not see this on the display, contact Orthofix Customer Service.) A backlight button is on the control unit. In low light, press the backlight button for illumination of the LCD.


Timing of Treatment Sessions

The Physio-Stim automatically times each treatment session. The timing begins when the device is turned on. The LCD shows a countdown of the time remaining in the treatment session. At the end of daily treatment, the device will turn itself off. To stop treatment prior to the end of a treatment session, simply press the On/Off button. To resume treatment, press the On/Off button again. The LCD will display the remaining treatment time.

Note: For the countdown to function correctly, treatment sessions should be greater than 60 minutes duration.

Charging/Recharging the Battery

The Physio-Stim is powered by a rechargeable lithium-ion battery pack. A power supply is provided with the device. Use only the Orthofix charging system to charge the battery. **Note:** The Physio-Stim battery will require charging prior to the first use.

To charge/recharge the battery, simply plug the barrel connector end of the power supply into the charger port located on the control unit. Plug the line cord securely into the power supply. Plug the line cord into any standard AC wall outlet. A fully discharged battery may require up to 4 hours to charge completely.

The Physio-Stim battery can be recharged at any time the device is not in use. **It is strongly recommended that the device be recharged after completing daily treatment.**


Note: The Physio-Stim will not deliver treatment while charging.

When the device is on, the Physio-Stim LCD will show a battery capacity symbol. A flashing battery outline, the symbol $b \frac{1}{2}$ and an audible beep indicate a battery low condition and that the battery needs to be charged. See "Visual and Audio Indicators" for more information.

Device Application


Models 3202 and 3303

Physio-Stim Models that are “U-shaped” are specifically designed for placement on a limb (e.g., tibia, femur, radius). These models may be worn over clothing, or over a cast or external fixation device if present.

To apply

1. Begin with the Physio-Stim unbuckled.
2. Place the Physio-Stim so that it is centered over the treatment site. It is suggested that the skin or cast (if present) be marked to facilitate placement of the device.
3. Bring the elastic strap around the limb and fasten the clasp.
4. If strap adjustment is needed, remove the device and tighten or loosen the strap until it feels secure and comfortable.

Model 3313

Physio-Stim Model 3313 may be placed on the collarbone (clavicle) or larger limb such as the thigh (femur). Model 3313 may be worn over clothing, or over a cast or external fixation device if present.

To apply on a collarbone


1. Fasten one end of the elastic strap to the end opposite the control unit.
2. Place the Physio-Stim so that it is centered over the treatment site (either left or right collarbone). The control unit should be in front and visible.
3. Bring the elastic strap around the body under the opposite arm and fasten to the end adjacent to the control unit.
4. If strap adjustment is needed, remove the device and tighten or loosen the strap until it feels secure and comfortable.

To apply on a limb

1. Use the short strap provided and attach to the device. Fasten one end of the elastic strap to either end of the device.
2. Place the Physio-Stim so that it is centered over the treatment site. Gently bend the device in a U-shape around the limb. It is suggested that the skin or cast (if present) be marked to facilitate placement of the device.
3. Bring the elastic strap around the limb and fasten the clasp.
4. If strap adjustment is needed, remove the device and tighten or loosen the strap until it feels secure and comfortable.

Models 3314L and 3314R

Physio-Stim Models 3314L or 3314R are intended for placement on either the left or right shoulder (proximal humerus).


To Apply

1. Fasten one end of the elastic strap to the end opposite the control unit.
2. Place the Physio-Stim so that it is centered over the treatment site (left or right shoulder). The control unit should be in front and visible.
3. Bring the elastic strap around the body under the opposite arm and fasten on the end of the device next to the control unit.
4. If strap adjustment is needed, loosen or tighten the strap until it feels secure and comfortable.

Model 3315


Physio-Stim Model 3315 is intended for placement on the hip (proximal femur).

To Apply

1. Fasten one end of the elastic strap to the device.
2. Place the Physio-Stim so that it is centered over the affected hip (left or right).
3. Bring the elastic strap around the body and fasten on the opposite end of the device.
4. If strap adjustment is needed, loosen or tighten the strap until it feels secure and comfortable.

Equipment Classification and Device Symbol Descriptions

Symbol	Meaning
	Attention - Refer to Instructions for Use
	Type BF Applied Part
	On/Off
	Backlight Button
	Storage Temperature Range
	Year of Manufacture for Active Device
	Charger Port

Equipment Classifications

- Internally powered equipment
- Type BF applied part
- IEC 529 enclosure rating: IPX0
- Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or nitrous oxide.
- Mode of operation: intermittent operation

The use of accessories other than those specified may result in increased emissions or decreased immunity of the device.

The battery charger is provided with a 3-wire appliance inlet but considered double insulated with Class II construction throughout.

For safe usage, follow manufacturer instructions when using the product.

Use of the product in any other manner could have harmful effects and/or void the warranty.

Note: Inspect the device prior to each use for wear or deterioration.
Do not use if the device does not appear to be in suitable condition.

Care and Cleaning

Physio-Stim is a technologically advanced electronic device and should be handled with appropriate care. Dropping or other mishandling of the Physio-Stim may cause damage to the device.

DO NOT expose the Physio-Stim to direct sunlight for long periods of time.

DO NOT expose the Physio-Stim to excessive heat. Avoid storing the device in areas prone to extreme temperatures such as an enclosed automobile or trunk.

DO NOT expose the Physio-Stim to excessive moisture.

DO NOT dispose of the Physio-Stim in an incinerator.

DO NOT use solvents to clean the Physio-Stim. Clean the device by wiping with a soft, damp cloth.

Travel

When traveling by air, it is best to check the Physio-Stim with the luggage. If the device is taken on board the airplane, it should not be worn when passing through passenger screening devices. The Physio-Stim could be damaged. The Physio-Stim user manual should be taken with you to quickly and easily identify the device for any security personnel.

Storage

Storage temperature range: -10°C to 45°C (14°F to 113°F)

Operating temperature range: $+5^{\circ}\text{C}$ to 40°C (41°F to 104°F)

Relative humidity: Up to 95%, noncondensing

Disposal

The Physio-Stim is for single patient use. Product contains lithium batteries; do not incinerate. Dispose of device properly to prevent injury. Please dispose of this product at collection facilities for waste electrical equipment used in household.

Service

If you have questions concerning the device or require any assistance, please call 800-535-4492 or 214-937-2000. There are no user serviceable parts. Notify the manufacturer for any servicing needs.

Warranty Policy

Orthofix Inc. warrants the Physio-Stim to be free from defects in materials and workmanship for one year from the date of first use. Provided that all terms and conditions of this Limited Warranty are complied with, Orthofix Inc. will replace defective components.

This Limited Warranty applies to the product only under normal use and does not cover any damage or defect caused by accident, misuse, abuse, fire, flood, and acts of God or by any alteration, tampering, repair or attempted repair by anyone other than Orthofix Inc. This warranty only applies to the patient for whom the product is prescribed and is not assignable or transferable. Defective products covered by this Limited Warranty must be returned to Orthofix Inc. Attention: Orthofix Returns. You must contact your local distributor to obtain the Return Authorization (RA) number and address prior to returning the product.

Except as specifically required by applicable law, the foregoing warranty is in lieu of all other warranties, expressed or implied and Orthofix Inc. specifically disclaims any and all warranties of merchantability or fitness for a particular purpose. Under no circumstances shall Orthofix Inc., its authorized representative, affiliated or subsidiary companies be liable for special, consequential or incidental damages. The sole remedy with respect to any defective product shall be limited to replacement.

This Limited Warranty may not be extended or modified except in writing by Orthofix Inc. No sales person, representative, distributor or physician is authorized to make or consent to any extension or modification of the terms of this Limited Warranty.


Manufactured by:
Orthofix Inc.
3451 Plano Parkway
Lewisville, TX 75056-9453
214-937-2000

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PS-1422 PL-US © Orthofix Holdings, Inc.