 **SPINAL-Stim[®]**
Osteogenesis Stimulator
U.S.A. EDITION

Table of Contents

Spinal-Stim® Instruction Manual

Prescription Information	1
Orthofix PEMF Stimulation	2
Clinical Success of the Spinal-Stim	2
Device Description	3
Treatment Instruction	4
Device Operation	4
Turning the Device On and Off	4
Timing of Treatment Sessions	5
Charging/Recharging the Battery	5
Visual and Audio Indicators	6
Device Application	7
Daily Application	9
Clinical Data Summary	11
Equipment Classification and Device Symbol Descriptions	12
Care and Cleaning	13
Travel	13
Storage	13
Disposal	13
Service	13
Warranty Policy	14

Package Contents:

- I- Spinal-Stim Osteogenesis Stimulator
- I- Literature Pack
- I- Power Supply
- I- Line Cord

THIS DEVICE IS NONSTERILE.
IT DOES NOT REQUIRE STERILIZATION.

U.S. Patent No. 5,743,844

Prescription Information

Indication

The Spinal-Stim® is a noninvasive electromagnetic bone growth stimulator indicated as a spinal fusion adjunct to increase the probability of fusion success and as a nonoperative treatment of salvage of failed spinal fusion, where a minimum of nine months has elapsed since the last surgery.

Contraindication

Cardiac pacemakers may be adversely affected by exposure to PEMF. Use of this device is contraindicated where the individual has an implanted cardiac pacemaker.

Warnings

- Although animal teratological studies performed with the device demonstrated no adverse findings, the safety of use of this device during pregnancy and nursing in humans has not been established.
- The safety and effectiveness of the use of this device on individuals lacking skeletal maturity have not been established.
- Animal studies conducted to date do not suggest any long-term adverse effects from the use of a similar device. However, long-term effects in humans are unknown.

Precautions

- This device should not be used if there are mental or physical conditions which preclude compliance with the physician and device instructions.
- This device has not been evaluated in treating patients with the following conditions: osseous or ligamentous spinal trauma, spondylitis, Paget's disease, moderate to severe osteoporosis, metastatic cancer, renal disease, and uncontrolled diabetes mellitus.
- The results of premarketing data from the randomized double-masked cohort indicate that inconsistent users (defined as those patients that used the device for less than an average of two hours per day) had success rates similar to those in the placebo group. Therefore, the use of the device for less than the minimum recommended usage may result in lower success rates.

Adverse Events

Rare instances of reversible minor discomfort have been reported. They were: cumbersome or uncomfortable, minor tingling or pain, minor skin rash, insomnia, fainting, nausea/diarrhea, and polymenorrhea.

Orthofix PEMF Stimulation

Pulsed electromagnetic field (PEMF) osteogenesis stimulation is a safe, nonsurgical, prescription treatment to promote spinal fusion and heal nonunion fractures. Electrical currents have been used to heal bones since the mid-1800s. However, it wasn't until the 1950s that scientists made an important discovery. When human bone is bent or broken, it generates an electrical field. This low-level electrical field activates the body's own repair mechanism which, in turn, stimulates bone healing.

Orthofix PEMF osteogenesis stimulators generate a uniform, low-level, pulsed electromagnetic field similar to the electrical field generated by the body. The application of PEMF directly to the fusion site helps activate and augment the body's natural healing process to enhance bone fusion. Thank you for including Orthofix in your healing process. To learn more about osteogenesis stimulation, please visit our website at www.bonestimulation.com.

Clinical Success of the Spinal-Stim


The Spinal-Stim has been proven safe and effective in clinical studies. A randomized, placebo-controlled clinical study demonstrated a 92% success rate (fusion rate) with the adjunctive use of Spinal-Stim. An open clinical trial of patients who were diagnosed with failed spinal fusions (nonoperative salvage) demonstrated that application of Spinal-Stim resulted in a 67% success rate, meaning 67% of the patient population in this study healed without further surgery. For additional information, see the "Clinical Data" section.

Device Description

The Spinal-Stim is an external, low-level, PEMF device and has been designed with patient comfort and convenience in mind. It is a single-piece device that is lightweight, flexible and portable, allowing freedom of movement during treatment. A Liquid Crystal Display (LCD) and audible alarm provide important feedback during treatment such as operational status, treatment time remaining, battery capacity, etc. See “Visual and Audio Indicators” for more information.

The Spinal-Stim is comprised of a control unit and treatment transducer. The control unit contains a micro-processor that generates the Spinal-Stim electrical signal. That signal is converted to a highly uniform, low-energy magnetic field by the treatment transducer. When the device is centered over the treatment area, the therapeutic Spinal-Stim PEMF signal is delivered directly to the treatment site.

Model: 2212

The Spinal-Stim is powered by a rechargeable lithium-ion battery pack. The LCD and audible alarm will alert the patient when the battery is low and needs to be recharged. See “Charging/Recharging the Battery” for more information. To ensure that the device is functioning properly, the Spinal-Stim constantly monitors battery voltage and the electrical signal. If at any time during treatment the device stops functioning properly, the LCD will display an appropriate symbol or error code. See “Visual and Audio Indicators” for more information.


Treatment Instructions

The Spinal-Stim should be worn daily for a minimum of two hours per day. Based upon clinical data, the overall treatment duration ranges between 90 and 270 days based upon specific patient conditions. At the end of the daily treatment prescribed by the physician, the device will turn itself off. The Spinal-Stim may be used at any time of day that is most convenient and comfortable for the patient. It is lightweight and adjustable. And because the Spinal-Stim is portable, treatment can be received while sitting, walking, reclining, sleeping, etc. However, since each patient is unique, the overall activity level should be based on physician instructions.

Device Operation

Turning the Device On and Off

The Spinal-Stim is turned on and off by pressing the On/Off button on the control unit of the device. When the device is on, a sequence of status messages will display momentarily. The LCD should then show the treatment time remaining and a flashing Orthofix logo. The flashing logo indicates that the device is on and functioning normally. (If you do not see this on the display, contact Orthofix Patient Services.) A backlight button is on the control unit. In low light, press the backlight button for illumination of the LCD.


Timing of Treatment Sessions

The Spinal-Stim automatically times each treatment session. The timing begins when the device is turned on. The LCD shows a countdown of the time remaining in the treatment session. At the end of daily treatment, the device will turn itself off. To stop treatment prior to the end of a treatment session, simply press the On/Off button. To resume treatment, press the On/Off button again. The LCD will display the remaining treatment time.

Note: For the countdown to function correctly, treatment sessions should be greater than 60 minutes duration.

Charging/Recharging the Battery

The Spinal-Stim is powered by a rechargeable lithium-ion battery pack. A power supply is provided with the device. Use only the Orthofix charging system to charge the battery. **Note:** The Spinal-Stim battery will require charging prior to the first use.

To charge/recharge the battery, simply plug the barrel connector end of the power supply into the charger port located on the control unit. Plug the line cord securely into the power supply. Plug the line cord into any standard AC wall outlet. A fully discharged battery may require up to 4 hours to charge completely.

The Spinal-Stim battery can be recharged at any time the device is not in use. It is strongly recommended that the device be recharged after completing daily treatment.

Note: The Spinal-Stim will not deliver treatment while charging.

When the device is on, the Spinal-Stim LCD will show a battery capacity symbol. A flashing battery outline, the symbol and an audible beep  indicate a battery low condition and that the battery needs to be charged. See "Visual and Audio Indicators" for more information.

Visual and Audio Indicators

The LCD and audible alarms are designed to provide helpful information to the user. The chart below shows the various displays and alarms and their meaning.

Spinal-Stim LCD Visual and Audio Indicators

Symbol / Alarm	Description	Meaning
	All LCD symbols visible and continuous audible alarm for approximately 5 seconds	Power-on self test
	Countdown timer displays remaining treatment time (hours & minutes) Orthofix logo flashes	Normal treatment in progress
	Countdown timer displays three dashes audible alarm (5 beeps)	No treatment time remaining
	Steady symbol for approximately 5 seconds	Treatment complete — power off
	Symbols flash / audible alarm (approximately 1 beep per second)	Battery low — recharge required
	Steady symbol indicates approximate % of charge Symbol filling repeatedly indicates charge mode	Battery status — remaining charge or charging mode
	Continuous audible alarm	Device locked — call for service
	Display of any E code (e.g., E01, E02 . . .)	Error message — call for service

Device Application


Putting on the Spinal-Stim for the First Time

Applying the Spinal-Stim for the first time is easy. And wearing the Spinal-Stim is comfortable if applied correctly. The Spinal-Stim may be worn over clothing, a brace or orthotic device. For first time use, follow the steps below.


Step 1:

Put the Spinal-Stim around the body to determine whether belt adjustment is needed. When properly adjusted, the device should be centered over the treatment site, and feel comfortable and secure. The control unit should be in front and on top as shown.


Step 2:

If significant belt adjustment is needed, lay the unit face down with the control unit on the table surface.


Step 3:

Release the latch.


Step 4:

Open the unit fully.


Step 5:


Open the Velcro® flap.


Step 6:
Pull the nonbuckled strap in or out to the desired length.

Step 7:
Store excess belt flat as shown.

Step 8:
Close the Velcro® flap and close the latch.


Step 9:
The length of the buckled strap and the nonbuckled strap should be approximately the same. If significant adjustment is needed, repeat the steps above.

Step 10:
For small adjustments, loosen or tighten the buckled strap. Your Spinal-Stim is now ready for daily application.


Daily Application

For daily application, clothing should be worn between the skin and the Spinal-Stim for greater comfort. The following is one suggested method for applying the device.


Step 1:

Begin with the Spinal-Stim unbuckled and with the control unit in front and on top.

Example of proper lumbar application

Step 2:

Hold both treatment transducers in one hand.


Step 3:

Grasp the back transducer by reaching behind with the opposite hand.



Step 4:

Fasten the buckle like a seatbelt. Make sure that the control unit is in front and on top.

Step 5:

Make any final adjustments by tightening or loosening the buckled strap.

Step 6:

The treatment transducers should be centered on the front and back of the body, over the fusion site.

Clinical Data Summary

Spinal-Stim has been tested in a clinical study involving 54 surgeons at 31 centers. This clinical investigation contained a prospective randomized double-masked trial of PEMF efficacy. Spinal-Stim was tested as a surgical adjunct in patients undergoing a first attempt at lumbar fusion. At one year postoperative, patients using active devices on a consistent daily regimen (an average of at least two hours per day) developed solid fusion in 92.2% of the cases.¹ Patients consistently using placebo (inactive) devices developed solid fusion in 67.9% of the cases. This 35% increase in treatment effect is statistically significant, and is realized regardless of:

- Number of Levels
- Gender
- Smoking
- Graft Type
- Vertebral Level
- Age
- Internal Fixation

The success rate for patients in the randomized double-masked phase for whom success or failure status is known at four years after treatment with the Spinal-Stim for all subjects (consistent and inconsistent users combined) was 63% (n=88) as compared with 83% in this phase of the clinical trial (i.e., one year postoperative). Consistent users (n=64) of the device in this phase had an initial success rate of 92.2% with a success rate of 70% after four years (a 24% reduction). Inconsistent users (n=34) and the entire placebo group (n=97) had an initial success rate of 65% with a success rate of 50% after four years (again, a 24% reduction). Long-term follow-up data indicates the success rate differentials between active and placebo units are maintained over time. Long-term, consistent Spinal-Stim users benefit with a 40% increase in fusion success, when compared to inconsistent and placebo device users. Based on this analysis, the reduction in long-term success rates appears unrelated to treatment with the Spinal-Stim. During this four year period, 10% of the original patients in the randomized double-masked phase were lost to follow-up and are not reflected in these success rates.


Spinal-Stim was also tested for nonoperative salvage in patients presenting with established pseudarthrosis of lumbar fusion in an open trial. Without concomitant regrafting, 67% of these cases reached a successful fusion with consistent (an average of at least 2 hours per day) PEMF treatment.² The four year success rates for these patients in the open trial, non-operative salvage phase for all subjects (consistent and inconsistent users combined) was 39% (n=119) as compared with 57% in this phase of the original clinical trial (i.e., one year postoperative). Consistent users (n=93) of the device in this phase had a success rate of 44% after four years. Inconsistent users (n=26) of the device in this phase had a success rate of 19% after four years. The reduction in success rates from the time of commercial marketing compared

with those at four years showed a similar percentage decrease (31%) to those in the randomized double-masked trial. During this four year period, 6% of the original patients in the open phase were lost to follow-up and are not reflected in these success rates.

1. Mooney, V., "A Randomized Double-Blind Prospective Study of the Efficacy of Pulsed Electromagnetic Field for Interbody Lumbar Fusions", SPINE, Vol. 15, No. 7, P708, 1990.

2. Simmons, JW, Hayes, MA, Christensen, KD, Dwyer, AP, Koulisis, CW, Kimmich, SJ: "The Effect of Postoperative Pulsed Electromagnetic Fields on Lumbar Fusion: Open Trial Study". Presented at the Annual Meeting of the North American Spine Society, Quebec City, Canada, 2 July 1989.

Equipment Classification and Device Symbol Descriptions

Symbol	Meaning
	Attention - Refer to Instructions for Use
	Type BF Applied Part
	On/Off
	Backlight Button
	Storage Temperature Range
	Year of Manufacture for Active Device
	Charger Port

Equipment Classifications

- Internally powered equipment
- Type BF applied part
- IEC 529 enclosure rating: IPX0
- Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or nitrous oxide.
- Mode of operation: intermittent operation

The use of accessories other than those specified may result in increased emissions or decreased immunity of the device.

The battery charger is provided with a 3-wire appliance inlet but considered double insulated with Class II construction throughout.

For safe usage, follow manufacturer instructions when using the product.

Use of the product in any other manner could have harmful effects and/or void the warranty.

Note: Inspect the device prior to each use for wear or deterioration.

Do not use if the device does not appear to be in suitable condition.

Care and Cleaning

Spinal-Stim is a technologically advanced electronic device and should be handled with appropriate care. Dropping or other mishandling of the Spinal-Stim may cause damage to the device.

DO NOT expose the Spinal-Stim to direct sunlight for long periods of time.

DO NOT expose the Spinal-Stim to excessive heat. Avoid storing the device in areas prone to extreme temperatures such as an enclosed automobile or trunk.

DO NOT expose the Spinal-Stim to excessive moisture.

DO NOT dispose of the Spinal-Stim in an incinerator.

DO NOT use solvents to clean the Spinal-Stim. Clean the device by wiping with a soft, damp cloth.

Travel

When traveling by air, it is best to check the Spinal-Stim with the luggage. If the device is taken on board the airplane, it should not be worn when passing through passenger screening devices. The Spinal-Stim could be damaged. The Spinal-Stim user manual should be taken with you to quickly and easily identify the device for any security personnel.

Storage

Storage temperature range: -10°C to 45°C (14°F to 113°F)

Operating temperature range: $+5^{\circ}\text{C}$ to 40°C (41°F to 104°F)

Relative humidity: Up to 95%, noncondensing

Disposal

The Spinal-Stim is for single patient use. Product contains lithium batteries; do not incinerate. Dispose of device properly to prevent injury. Please dispose of this product at collection facilities for waste electrical equipment used in household.

Service

If you have questions concerning the device or require any assistance please call 800-535-4492 or 214-937-2000. There are no user serviceable parts. Notify the manufacturer for servicing needs.

Warranty Policy

Orthofix Inc. warrants the Spinal-Stim to be free from defects in materials and workmanship for one year from the date of first use. Provided that all terms and conditions of this Limited Warranty are complied with, Orthofix Inc. will replace defective components.

This Limited Warranty applies to the product only under normal use and does not cover any damage or defect caused by accident, misuse, abuse, fire, flood, and acts of God or by any alteration, tampering, repair or attempted repair by anyone other than Orthofix Inc. This warranty only applies to the patient for whom the product is prescribed and is not assignable or transferable. Defective products covered by this Limited Warranty must be returned to Orthofix Inc. Attention: Orthofix Returns. You must call a Patient Services Representative at 1-800-535-4492 or your local distributor to obtain the Return Authorization (RA) number and address prior to returning the product.

Except as specifically required by applicable law, the foregoing warranty is in lieu of all other warranties, expressed or implied and Orthofix Inc. specifically disclaims any and all warranties of merchantability or fitness for a particular purpose. Under no circumstances shall Orthofix Inc., its authorized representative, affiliated or subsidiary companies be liable for special, consequential or incidental damages. The sole remedy with respect to any defective product shall be limited to replacement.

This Limited Warranty may not be extended or modified except in writing by Orthofix Inc. No sales person, representative, distributor or physician is authorized to make or consent to any extension or modification of the terms of this Limited Warranty.

Manufactured by:
Orthofix Inc.
345 I Plano Parkway
Lewisville, TX 75056-9453
214-937-2000

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SS-1409 PL-US © Orthofix Holdings, Inc.